


Stop Bleeding Fast with the Quick-Tight Tourniquet



91160
Quick-Tight Tourniquet

- Quick tightening with a single strap attachment twist point on the windlass wand
- Strong solid windlass wand with etched ends for performance grip
- Dual method of applying over wound with detachable D-Clip. User can unhook clasp and wrap around limb above wound or keep strap clasped and slide up limb above wound (see diagram below)
- Tourniquet won't loosen with proven hook and loop backing on the strap
- Wide 1.5" constricting band to increase constricting pressure and provide patient comfort
- "Time" label informs medical personnel when tourniquet was applied

Dual Methods of Applying Over a Wound



Unhook the clasp, loosen strap, wrap around limb, reattach clasp.

OR Slide strap above wound.